

WBtool

Create and annotate composite Western blot figures from pre-rotated gel images

1. Install the Script

1. Open link and install Fiji: <https://imagej.net/software/fiji/downloads>

2. Input WBtool plugin into Fiji

a. Open link and download wbtool_fiji_2.0.py: <https://github.com/masha-rgfj/wbtool-fiji-plugin>

b. Install plugin, two ways:

i. Fiji file > Plugins > drag wbtool file into Plugins folder

ii. Open the Fiji app. Plugins > Install > open wbtool file > save to plugins

2. Use the Tool

1. In Fiji, select Plugins > wbtool

2. Basics: Upload Image/s, Mark kDa, & Crop (must do in order)

a. , upload image from device. It will open in a new window.

b. , click on band, type in value, repeat for all bands

c. , drag mouse over image to crop region. Name protein.

i. to add another image

ii. and to adjust size

iii. and to adjust placement (multiple images)

3. Label samples:

a. , click on band, add label and angle, repeat for all lanes

4. Special Annotations:

a. add horizontal line

b. label specific bands

c. add text anywhere

5. Edit Annotations:

a. edit text

b. delete

c. change text size

d. change placement

6. Export: or

contact us

Masha Pirozhkova, WBtool creator
pirozhkova@mail.tau.ac.il
LCRB, Tel Aviv University

try it yourself!

